

Frusto-conical bottom Vessel Mini-Bioreactor

Mini-Bioreactor

For the best oxygen transfer

- ◆ The material used to form the frusto-conical bottom vessel mini-bioreactor was specially chosen for its ideal physical and chemical surface nature (see SEM micrograph), which provides optimal dissolved oxygen generation as the culture medium sweeps upon the surface (Bioprocess International, June 2008).
- ◆ The frusto-conical bottom of the culture vessel on an orbital shaker causes the culture medium surface to form a tilted surface that sweeps the vessel's other, air-exposed surface (Bioprocess International, June 2008; WO2007/142664). This action generates dissolved oxygen faster than that of any of sparging or bubbling method and has the added advantage of eliminating bubble burst-caused cell damage (WO2007/142664). The frusto-conical bottom mini-bioreactor series (45ml, 150ml, 1000ml and 2000ml) is ideal for the process development and scale up seed train support of mammalian cells, insect cells, plant cells, yeast cells and bacterial cells.
- ◆ Successful cell cultures using the above mini-scale vessels are supported by scale-up bioreactors (5 liter, 50 liter, 150 liter and 300 liter) both for suspension and paper-disc carrier perfusion cultures (see www.amprotein.com or www.amprotein-china.com).



Vessel surface (micrograph)

Frusto-conical Vessel



| Size | Cell Type | Culture Method | Highest Cell Density | Media | Optimization |
|-------|----------------|----------------|-------------------------------------|----------|--------------|
| 40ml | CHO-S | Fed-batch | $1.8 \sim 2.5 \times 10^7$ cells/ml | B001 | None |
| | CHO-K1 | Fed-batch | $1.2 \sim 1.5 \times 10^7$ cells/ml | B001 | None |
| | SF9 (insect) | Batch | 1.8×10^7 cells/ml | SFM90 II | None |
| | Tobacco Cell | Batch | PCV = 75~80% | MS | Optimized |
| | Carrot Callus | Batch | ↑ Growth 20 - 25 times | MS | Optimized |
| 150ml | E. coli | Batch | $OD_{600} = 2.2 \sim 2.5$ | LB | None |
| | CHO-S | Fed-batch | $1.5 \sim 2.0 \times 10^7$ cells/ml | B001 | None |
| 1L | CHO-K1 | Fed-batch | $1.2 \sim 1.4 \times 10^7$ cells/ml | B001 | None |
| | CHO-S | Fed-batch | $1.3 \sim 1.7 \times 10^7$ cells/ml | B001 | None |
| 2L | CHO-K1 | Fed-batch | $0.8 \sim 1.3 \times 10^7$ cells/ml | B001 | None |
| | CHO-S | Fed-batch | $1.3 \sim 1.5 \times 10^7$ cells/ml | B001 | None |
| | CHO-K1 | Fed-batch | $0.8 \sim 1.2 \times 10^7$ cells/ml | B001 | None |



Hangzhou AmProtein Bioengineering Co.,Ltd

Add : No.63 Jiuhuan Rd., Jianggan Technological

& Economic Park, Hangzhou, Zhejiang, China 310019

Tel : 86-571-28926718/28926719 Fax : 86-571-28925102

Web : www.amprotein.com / www.amprotein-china.com

Email : info@amprotein.com